



Seeplex®

Detection of Respiratory Pathogens

Influenza A/B virus typing
7 Respiratory virus detection
12 Respiratory virus detection
15 Respiratory virus detection
6 Respiratory bacteria detection

CE IVD
www.seegene.com

 **Seegene**



Seeplex[®] Respiratory Pathogen Detection Products

Specimen

Nasopharyngeal swab
Nasopharyngeal aspirate
Bronchoalveolar lavage
Sputum (for detection of respiratory bacteria)

Influenza A/B OneStep Typing

Seeplex[®]

ACE

Medical
Device
License by

Health
Canada

CE

IVD

Features

- Typing of influenza A and B viruses
- Subtyping of influenza A (H1, H3)
- Differentiation of (H1N1) & seasonal influenza A (H1N1)
- cDNA synthesis and multiplex PCR in one-step reaction

WHO information for laboratory diagnosis of
pandemic (H1N1) 2009 virus in humans (18 August 2009)

Influenza A virus (Common)
2009 pandemic H1
Human seasonal influenza A (H1)
Human seasonal influenza A (H3)

Follow
WHO
Guideline

Analytes - Set (A)


Influenza B virus
Influenza A virus (Common)
2009 pandemic H1
Human seasonal influenza A (H1)
Human seasonal influenza A (H3)

Interpretation of laboratory results :


PCR -A sample is considered positive if results from tests using two different PCR targets (e.g. primers specific for universal M gene and swine H1 haema- gglutinin gene) are positive but the PCR for human H1 + H3 is negative. If RT-PCR for multiple haemagglutinin (HA) targets (i.e. H1, H3, and H1-pandemic) give positive results in the same specimen, the possibility of PCR contamination should first be excluded by repeating PCR procedure using new RNA extract from the original specimen or RNA extract from another specimen. If repeated positive results for multiple HA targets are obtained, this raises the possibility of co-infection, which should be confirmed by sequencing or virus culture.

(http://www.who.int/csr/resources/publications/swineflu/WHO_Diagnostic_Recommendations_H1N1_20090521.pdf)

* 96 well format



• Seeplex[®] Influenza A/B OneStep Typing (FluB, FluA, 2009 pandemic H1, H1, H3)



94 samples
470 results
in 5.5h

• Existing Real-time RT-PCR (FluA, 2009 pandemic H1, H1, H3)



44 samples
176 results
in 5h

RV12 ACE Detection

Simultaneous detection of 6 respiratory bacteria

Seeplex®

ACE

Medical
Device
License by

Health
Canada

CE

IVD

Result | ScreenTape® System (Agilent Technologies)

1~6: Clinical samples, N: Negative control

Set (A)

Lane	L2	L3	L4	L5	L6	L7	L8
Sample ID	1	2	3	4	5	6	N
Internal Control	+	98	+	101	+	99	+
MPV	+	96	+	92	+	94	+
AdV	+	96	+	92	+	94	+
229E/NL63	+	96	+	92	+	94	+
RSV A	+	20	+	20	+	57	+
RSV B	+	20	+	20	+	57	+
Flu A	+	20	+	20	+	57	+
Flu B	+	20	+	20	+	57	+
Unidentified	+	20	+	20	+	57	+

Set (B)

Lane	L2	L3	L4	L5	L6	L7	L8
Sample ID	1	2	3	4	5	6	N
Internal Control	+	98	+	96	+	98	+
Flu B	+	98	+	96	+	98	+
OC43	+	98	+	96	+	98	+
HRV A/B	+	98	+	96	+	98	+
RSV A	+	98	+	96	+	98	+
Flu A	+	98	+	96	+	98	+
RSV B	+	98	+	96	+	98	+
Unidentified	+	98	+	96	+	98	+

Analytes - Set (A)

Metapneumovirus
Adenovirus A/B/C/D/E
Coronavirus 229E/NL63
Parainfluenza virus 1,2,3

Analytes - Set (B)

Influenza virus A/B
Respiratory syncytial virus A, B
Rhinovirus A/B
Coronavirus OC43

RV7 Detection

Simultaneously detecting major respiratory viruses

Seeplex®

GEL

Medical
Device
License by

Health
Canada

CE

IVD

A | Agarose gel

Lane	1	2	3	4	5	6	7	8	9	10	N
Sample ID	1	2	3	4	5	6	7	8	9	10	N
Internal Control	+	+	+	+	+	+	+	+	+	+	+
Flu B	+	+	+	+	+	+	+	+	+	+	+
AdV	+	+	+	+	+	+	+	+	+	+	+
MPV	+	+	+	+	+	+	+	+	+	+	+
RSV	+	+	+	+	+	+	+	+	+	+	+
Rhino A	+	+	+	+	+	+	+	+	+	+	+
P1V	+	+	+	+	+	+	+	+	+	+	+
Flu A	+	+	+	+	+	+	+	+	+	+	+
Unidentified	+	+	+	+	+	+	+	+	+	+	+

Analytes

Influenza B virus
Adenovirus
Metapneumovirus
Respiratory syncytial virus A/B
Rhinovirus A
Parainfluenza virus 1/2/3
Influenza A virus

Related publications

- Kim et al., Rapid detection and identification of 12 respiratory viruses using a dual priming oligonucleotide system-based multiplex PCR assay. J Virol Methods. 2009 Mar;156(1-2): 111-6. Epub 2008 Dec 247.
- Dreus et al., Use of the Seeplex RV Detection kit for surveillance of respiratory viral outbreaks in Toronto, Ontario, Canada. Ann Clin Lab Sci. 2008;38(4):376-379.
- Kim et al., Detection and clinical manifestations of twelve respiratory viruses in hospitalized children with acute lower respiratory tract infections : Focus on human metapneumovirus, human rhinovirus and human coronavirus. Korean Journal of Pediatrics, 2008; 51(8):834-841.

PneumoBacter ACE Detection

Simultaneous detection of 6 respiratory bacteria

Seeplex®

ACE

CE

IVD

Result | ScreenTape® System (Agilent Technologies)

1~7: Clinical samples, N: Negative control

Lane	L2	L3	L4	L5	L6	L7	L8
Sample ID	1	2	3	4	5	6	N
Internal Control	+	97	+	96	+	97	+
M.pneumoniae	+	48	+	48	+	48	+
L.pneumoniae	+	48	+	48	+	48	+
S.pneumoniae	+	48	+	48	+	48	+
H.influenzae	+	26	+	26	+	26	+
B.pertussis	+	26	+	26	+	26	+
C.pneumoniae	+	90	+	90	+	90	+
Unidentified	+	90	+	90	+	90	+

Bacteria causing typical pneumonia

Streptococcus pneumoniae
Haemophilus influenzae

Bacteria causing atypical pneumonia

Chlamydomphila pneumoniae
Legionella pneumophila
Bordetella pertussis
Mycoplasma pneumoniae

Related publications

- Cho MC, et al., Ann Lab Med, 2012;32: 133-138
- Loes K, et al., J Microbiol Methods, 2010; 82(2): 131-5.
- Lim JS et al., J Microbiol Methods, 2010; 53(3): 373-9.
- Jung CL, et al., Korean J Clin Microbiol, 2010;13(1): 40-6.
- Shin JH, et al., Korean J Pediatr, 2009; 52(12): 1358-63.
- Lee KS, et al., Korean J Pediatr Infect Dis, 2009;16 (2):175-82.
- Park J, et al., Korean J Lab Med, 2009; 29(4): 307-13.
- Higgins RR, et al., Ann Clin Microbiol Antimicrobial, 2009;8(10):1-9.

Seeplex® Automatic System



Automated Extraction and PCR Setup

Nimbus IVD

Certification



Automated Liquid Handling Workstations from NA Extraction to PCR Setup

- | | |
|--|--|
| <ul style="list-style-type: none"> A. Automated DNA/RNA extraction and PCR setup B. Wide range of pipetting volume :
0.5 ~ 1,000 uL C. Minimized hands-on time D. No manual pipetting required | <ul style="list-style-type: none"> E. 48 samples processed in 2.5hrs F. Enclosed system and liquid waste system preventing contamination G. Built in UV lamp for decontamination of residual DNA/RNA H. Handheld barcode scanner |
|--|--|



(Hamilton)

Cat. No. 65415-02

Automatic Detection & Analysis

LabChip® Dx

Certification



Accelerating throughput on your PCR work

- | | |
|---|--|
| <ul style="list-style-type: none"> A. Convenient and no post-PCR contamination <ul style="list-style-type: none"> - No need for sample dilution and transfer (ex. Bio-RAD Hard-shell 96 well plate) B. High throughput <ul style="list-style-type: none"> - Minimum 1hr for 96-well plate C. Precise and reliable data
Stable reagents | <ul style="list-style-type: none"> D. Improvements on time and labor savings <ul style="list-style-type: none"> - Powerful automated electrophoresis system with barcode reader E. Barcode system for effective sample traceability F. Comprehensive and user friendly software |
|---|--|




(Caliper LifeSciences)


Cat. No. 122000

Sample Collection Method

Nasopharyngeal aspirates



Nasopharyngeal swab



- **Nasopharyngeal aspirates**

Prepare a bulb syringe or suction apparatus and 3~5ml of saline to collect the specimen. Squeeze the bulb to inject saline, then release the bulb to suction it out.

- **Nasopharyngeal swab**

Stir the brush into the bottle containing 1~3 ml buffer solution (saline or 1X PBS).

- **Bronchoalveolar lavage (BAL)**

Method used for a cleaning purpose during the bronchoscopy process.

RV15 ACE Detection

Simultaneous detection of currently the most known respiratory viruses

Seeplex®

ACE

Medical
Device
License by
Health
Canada

CE

IVD

Analytes - Set (A)

Adenovirus A/B/C/D/E
Coronavirus 229E/NL63
Parainfluenza virus 1
Parainfluenza virus 2
Parainfluenza virus 3

Analytes - Set (B)




Coronavirus OC43
Rhinovirus A/B/C
Influenza A virus
Respiratory syncytial virus A
Respiratory syncytial virus B

Analytes - Set (C)

Bocavirus 1/2/3/4
Influenza B virus
Metapneumovirus
Parainfluenza virus 4
Enterovirus

Result | ScreenTape® System (Agilent Technologies)

1~5 : Clinical samples, N : Negative control



RV15 OneStep ACE Detection

Seeplex®

ACE

Medical
Device
License by
Health
Canada

CE

IVD

- Including whole process control(extraction & RT-PCR control)
- cDNA synthesis and multiplex PCR in one-step reaction

Analytes - Set (A)

PCR Control
Adenovirus A/B/C/D/E
Coronavirus 229E/NL63
Parainfluenza virus 1
Parainfluenza virus 2
Parainfluenza virus 3

Analytes - Set (B)




PCR Control
Coronavirus OC43
Rhinovirus A/B/C
Influenza A virus
Respiratory syncytial virus A
Respiratory syncytial virus B

Analytes - Set (C)

Whole Process Control
Bocavirus 1/2/3/4
Influenza B virus
Metapneumovirus
Parainfluenza virus 4
Enterovirus

Result | ScreenTape® System (Agilent Technologies)

1~7 : Clinical samples, N : Negative control



⚠ : Only for set C negative control, no band should appear for WPC.

Seeplex® Key Products

Infectious Pathogen Detection

1. Respiratory Pathogen Detection
2. Tuberculosis Detection
3. STD Pathogen Detection
4. Human Papillomavirus Detection
5. HSV2 ACE Detection
6. Meningitis ACE Detection
7. Diarrhea ACE Detection

Somatic Mutation Detection

1. BRAF V600E ACE Detection
2. Leukemia Detection

Drug Resistance Detection

1. ClaR-H. pylori ACE Detection
2. HBV Lami-resistant ACE
3. VRE ACE Detection

SNP Detection

1. ApoE Genotyping
2. CYP2C19 ACE Genotyping
3. MTHFR Genotyping

Respiratory Pathogen Detection Products

	Cat. No.	Size
Seeplex® RV7 Detection	RV3210	50 rxns
Seeplex® RV12 ACE Detection	RV6C00Y	50 rxns
Seeplex® RV15 ACE Detection	RV6F00Y	50 rxns
Seeplex® RV15 OneStep ACE Detection	RV6F01Y	50 rxns
Seeplex® Influenza A/B OneStep Typing	RV6520X	100 rxns
Seeplex® PneumoBacter ACE Detection	PB1610Y	50 rxns



Seegene, Inc.

Taewon Bldg., 65-5, Bangyi-Dong, Songpa-Gu, Seoul, Korea
TEL : +82-2-2240-4000, FAX : +82-2-2240-4046
E-mail: korea@seegene.com

Seegene USA

9801 Washingtonian Blvd. Suite 290, Gaithersburg, MD 20878, U.S.A.
TEL : +1-301-762-9066, FAX : +1-301-762-9088
E-mail: usa@seegene.com

Seegene Germany

Ludwig-Erhard Str. 30-34, D 65760, Eschborn, Germany
TEL : +49 (0) 6196-769-3468, FAX : +49 (0) 6196-769-4137
E-mail: eu@seegene.com

Not Available in U.S.